

# Mesothelin Control

# 360-20

This Lung Marker Control is an assayed, two-level control containing Cyfra and ProGRP antigens. The Fujirebio Diagnostics Lung Marker Control can be utilized as an independent, internal quality control to monitor the precision of the laboratory's Lumipulse G Cyfra and ProGRP assays. Laboratories can feel confident releasing their patient results when using a control to assess the long-term accuracy and precision of their Cyfra and ProGRP assays.

## Характеристики

**Control:** 1 vials of level 1  
1 vials of level 2

**Volume:** 2x1 mL

## Информация для заказа

Кат № 360-20  
Mesothelin Control

