

**Открытый ротор, 4-местный,
для пробирок по ASTM и склянок Бабкока (жиромер) 50 мл**



± 90°

$n = 2,000 \text{ мин}^{-1}$
макс. RCF 984

Кат. No. (без подвесов) **4474**

объем, мл	50	100	100	100	50
Ø x L, мм	45x130	58x161	44x168	37x200	36.5x185
Кат. No.	-	0528	-	0531	-
Кат. No. 4275  ротор Кат. No. 4474					
					
Cat. No.	4278-A	0771	4277	4276	0703
Ø отверстия x L, мм	-				
кол. пробирок на ротор	4				
макс. RCF	984	961		912	
радиус, мм	220	215		204	
время разгона, сек	17				
время остановки с торможением, сек	20				
температура °C	-8				