

Lung Marker Control

240-20

This Lung Marker Control is an assayed, two-level control containing Cyfra and ProGRP antigens. The Fujirebio Diagnostics Lung Marker Control can be utilized as an independent, internal quality control to monitor the precision of the laboratory's Lumipulse G Cyfra and ProGRP assays. Laboratories can feel confident releasing their patient results when using a control to assess the long-term accuracy and precision of their Cyfra and ProGRP assays.

Характеристики

Control: 2 vials of level 1
2 vials of level 2

Volume: 4x3 mL

Информация для заказа

Кат № 240-20
Lung Marker Control

