

**Угловой ротор, 14-местный**




α 36°

n = 9,000 / 10,000 мин<sup>-1</sup>  
макс. RCF 14,127 17,441

с биогерметизацией

**Кат. No. 5655**

объем, мл	15	25	7.5–8.5	9–10	10	8.5–10	15	50	30
Ø x L, мм	17 x 100	24 x 100	15 x 92	16 x 92	15 x 102	16 x 100	17 x 120	29 x 115	26 x 95
<b>Кат. No.</b>	<b>0518</b>	<b>0519</b>	<b>пробирки для крови/мочи</b>				<b>0509</b>	<b>0513</b>	<b>0545</b>
крышка с биогерметизацией вкл. 									
<b>ротор</b> Кат. No. <b>5655</b>									
<b>Кат. No.</b>	<b>1473</b>	<b>1472</b>	<b>1473</b>				<b>1474</b>	-	<b>1472</b>
Ø отверстия x L, мм	17.5 x 91.5	26 x 85	17.5 x 91.5				17 x 107	30 x 100	26 x 85
кол. пробирок на ротор	14								
макс. RCF классическая центрифуга/ с охлаждением	13,493 / 16,658	13,131 / 16,211	13,493 / 16,658				13,765 / 16,994	14,127 / 17,441	13,131 / 16,211
радиус, мм	149	145	149				152	156	145
время разгона, сек классическая центрифуга/ с охлаждением	124 / 145								
время остановки с торможением, сек классическая центрифуга/ с охлаждением	94 / 108								
температура, °C	+16								