

Угловой ротор, 30-местный




± 45°

n = 15,000 мин⁻¹
макс. RCF 24,400

с биогерметизацией,
фенолоустойчивый

Кат. No. 4489-A

объем, мл	0.2	0.4	0.5	0.8	1.5	2.0	0.5
Ø x L, мм	6x18	6x45	8x30	8x45	11 x38		10.7 x36
Кат. No.	-				2078	0536	Pediatric
крышка с биогерметизацией вкл.  ротор Кат. No. 4489-A							
							
Кат. No.	2024		2023		2031	-	0788
Ø отверстия x L, мм	6x40		8x40		10.2x19	11.2 x40.1	11.2 x40.1
кол. пробирок на ротор	30						15
макс. RCF	24,400						23,394
радиус, мм	97						93
время разгона, сек	73						
время остановки с торможением, сек					77		
температура °C	+7						