

CanAg ProGRP EIA

220-10

ProGRP is the precursor of the gut hormone GRP (Gastrin Releasing Peptide). Fragments of ProGRP in serum has been shown to be elevated in patients diagnosed with Small Cell Lung Cancer (SCLC). The assay is a one-step immunoassay that utilizes biotinylated Anti-ProGRP polyclonal antibody and horseradish peroxidase (HRP) labeled Anti-ProGRP monoclonal antibody E146 in Streptavidin coated microtiter strips for the measurement of ProGRP in serum.

SPECIFICATIONS

- Results within:** 3 hours, one step procedure
- Detection limit:** < 10 ng/L
- Measuring range:** 10-2 000 ng/L
(may be extended by sample dilution)
- Sample volume:** 50 µL
- Hook effect:** No hook up to 1 500 000 ng/L
- Incubation temp:** 20-25° C
- Detection:** 450 nm

ORDERING INFORMATION

- Prod. No. 220-10
- CanAg ProGRP EIA
- For 96 determinations*

