

Откидной ротор, 4-местный

± 90°

n = 4,000 мин⁻¹
max. RCF 3,095

изображен с подвесами 5051
и крышками 5053

Кат. No. (без подвесов) **5694**

объем, мл	1.5	2.0	1.5	2.0	5	6	7	9	15	25	50	100	100	1.1–1.4	2.6–3.4	2.7–3	4–5.5
Ø x L, мм	11 x 38		11 x 38		12 x 75	12 x 82	12 x 100	14 x 100	17 x 100	24 x 100	34 x 100	44 x 100	40 x 115	8 x 66	13 x 65	11 x 66	15 x 75
Кат. No.	-	-	2078	0536	0553	0501	0578	0500	0518	0519	0521	0526	0523	емкости для забора крови			
крышка Кат. No. 5053  подвес Кат. No. 5051																	
Кат. No.	5257		5281		5227		5247	5264	5248	5242	5243	5262	5249	5267	5268	5227	5264
Ø отверстия x L, мм	12.5x42		11.5x50		12.5x42		12.5x86	16x50	17.5x86	26x86	36x86	45.5x86	42x86	11 x 44	13.5x52	12.5x42	16x50
кол. пробирок на ротор	160		64		80		80	48	48	20	8	4		80	48	80	48
макс. RCF верх. / нижн.	1,950 / 2,826		2,826		2,773		2,755	2,773	2,755			2,755		2,737	2,808	2,773	
радиус, мм верх. / нижн.	109 / 158		158		155		154	155	154			154		153	157	155	
время разгона, сек	42																
время остановки с торможением, сек	51																
температура °C	-3																

объем, мл	4.5–5	4.9	7.5–8.5	9–10	10	1.6–5	4–7	4–7	8.5–10	15	50	12	30	50	двухуровневый адаптер для цитологических принадлежностей	
Ø x L, мм	11 x92	13x90	15x92	16x92	15x102	13x75	16x75	13x100	16x100	17x120	29x115	17x100	25x110	29x115		
Кат. No.	емкости для забора крови / пробирки для мочи										0509	0513	-	-	цитологические камеры	
крышка Кат. No. 5053  подвес Кат. No. 5051																
																
Кат. No.	5227	5268	5264	5258		5268	5264	6301	5248	6306	5259	6306	5266	5243	5280	
Ø отверстия x L, мм	12.5x42	13.5x52	16x50	17.6x86		13.5x52	16x50	13.5x86	17.5x86	17x90	30x86	17x90	26x86	36x86	-	
кол. пробирок на ротор	80	48	48	44		48				28	8	28	20	8	8	
макс. RCF верх. / нижн.	2,773	2,808	2,773	2,755		2,808	2,773	2,755		2,898	2,826	2,898	2,755		1,735/ 2,737	
радиус, мм верх. / нижн.	155	157	155	154		157	155	154		162	158	162	154		97/ 153	
время разгона, сек	42															
время остановки с торможением, сек	51															
температура °С	-3															



изображен с подвесами 5092

Откидной ротор 5694


объем, мл	5	6	7	10	15	25	50	100	100	250	1.1–1.4	1.1–1.4	2.6–3.4
Ø x L, мм	12x75	12x82	12x100	17x70	17x100	24x100	34x100	44x100	40x115	65x115	8x66	8x66	13x65
Кат. No.	0553	0501	0578	2079	0518	0519	0521	0526	0523	0530	емкости для забора крови		
крышка с биогерметизацией Кат. No. 5093  подвес Кат. No. 5092 													
													
Кат. No.	5128		5120	5136		5122	5124	5125	5126	1791	5137	5138	5137
Ø отверстия x L, мм	12.8x45		12.8x82	17.5x60.7		25.5x82	35.5x82.5	45.5x85	42x85	66x103	13.5x60.7	12.8x54.5	13.5x60.7
кол. пробирок на ротор	48			32		16	4				32	48	32
макс. RCF	3,005			2,952		2,898	2,952	2,952		3,095	2,952	2,540	2,952
радиус, мм	168			165		162	165	165		173	165	142	165
время разгона, сек	42												
время остановки с торможением, сек	51												
температура °C	-3												

объем, мл	2.6–3.4	2.7–3	2.7–3	4–5.5	4.5–5	4.9	7.5–8.5	9–10	10	1.6–5	1.6–5	4–7
Ø x L, мм	13x65	11 x66	11 x66	15x75	11 x92	13x90	15x92	16x92	15x102	13x75	13x75	13x100
Кат. No.	емкости для забора крови / пробирки для мочи											
крышка с биогерметизацией Кат. No. 5093  подвес Кат. No. 5092 												
												
Cat. No.	5138	5137	5138	5136	5120	5137	5136			5137	5138	5120
Ø отверстия x L, мм	12.8x54.5	13.5x60.7	12.8x54.5	17.5x60.7	12.8x82	13.5x60.7	17.5x60.7			13.5x60.7	12.8x54.5	12.8x82
кол. пробирок на ротор	48	32	48	32	48	32				48		
макс. RCF	2,540	2,952	2,540	2,952	3,005	2,952				2,540	3,005	
радиус, мм	142	165	142	165	168	165				142		168
время разгона, сек	42											
время остановки с торможением, сек	51											
температура °C	-3											



изображен с подвесами 5092

Откидной ротор 5694

объем, мл	4—7	8	8.5—10	15	50	12	25	30	50	10	250	290
Ø x L, мм	16x75	16 x 125	16x100	17x120	29x115	17 x 100	25x90	25x110	29x115	16x80	61x122	62x137
Кат. No.	емкости для забора крови			0509	0513	-					5127	-
крышка с биогерметизацией Кат. No. 5093  подвес Кат. No. 5092												
												
Кат. No.	5136	5121	5136	5129	5123	5121	5134	5122	5135	5136	6319	
Ø отверстия x L, мм	17.5x60.7	17.5x82	17.5x60.7	17x85	30x85	17.5x82	26x73	25.5x82	30x99	17.5x60.7	66x103	
кол. пробирок на ротор	32	28	32	28	8	28	12	16	8	32	4	
макс. RCF	2,952	3,059	2,952	3,095		3,005	2,826	2,898	3,023	2,952	3,095	
радиус, мм	165	171	165	173		168	158	162	169	165	173	
время разгона, сек	42											
время остановки с торможением, сек	51											
температура °C	-3											