

Лоток для инструментов KU 312 / KU 313, 200 100 50 мм, с ручкой / без ручки

ИМ-7.107461 до 22.04.2024

сталь, прямоугольный, для инструментов, длина 200 мм, с крышкой

Nora instruments Medizintechnik GmbH, Германия



| Länge | | |
|-----------------------|----------|--|
| 18,5 x 8,0 x 3,5 cm | KU 31216 | |
| 20,0 x 10,0 x 5,0 cm | KU 31220 | |
| 22,0 x 15,0 x 5,0 cm | KU 31222 | |
| 24,0 x 15,0 x 5,0 cm | KU 31226 | |
| 26,5 x 20,5 x 5,0 cm | KU 31230 | |
| 30,0 x 15,0 x 5,0 cm | KU 31231 | |
| 30,0 x 20,0 x 10,0 cm | KU 31250 | |
| 30,0 x 20,0 x 12,0 cm | KU 31251 | |

| Dicke der Klinge | | |
|------------------|----------|--|
| 1,0 mm | KU 31216 | |
| 1,5 mm | KU 31220 | |
| 2,0 mm | KU 31222 | |
| 2,5 mm | KU 31226 | |
| 3,0 mm | KU 31230 | |
| 3,5 mm | KU 31231 | |
| 4,0 mm | KU 31250 | |
| 4,5 mm | KU 31251 | |