

Лоток для инструментов KU 312 / KU 313, 305 205 50 мм, с ручкой / без ручки

ИМ-7.107461 до 22.04.2024

сталь, прямоугольный, для инструментов, длина 305 мм, с крышкой

Nora instruments Medizintechnik GmbH, Германия



Länge		
16,5 x 8,0 x 3,5 cm	KU 31216	
20,0 x 10,0 x 5,0 cm	KU 31220	
22,0 x 15,0 x 5,0 cm	KU 31222	
16,5 x 15,0 x 5,0 cm	KU 31216B	
20,5 x 20,5 x 5,0 cm	KU 31230	
10,0 x 15,0 x 5,0 cm	KU 31231	
10,0 x 20,0 x 10,0 cm	KU 31250	
10,0 x 20,0 x 12,0 cm	KU 31251	

Deckel ohne Klinge!		
Länge ohne Klinge		
Tiefe mit Griff		
Einseitig oder beidseitig		
Länge ohne Griff		
16,5 x 8,0 x 3,5 cm	KU 31316	
20,0 x 10,0 x 5,0 cm	KU 31320	
22,0 x 15,0 x 5,0 cm	KU 31322	
20,0 x 15,0 x 5,0 cm	KU 31326	
10,5 x 20,5 x 5,0 cm	KU 31330	
10,0 x 15,0 x 5,0 cm	KU 31331	
10,0 x 20,0 x 10,0 cm	KU 31350	
10,0 x 20,0 x 12,0 cm	KU 31351	