

Лоток для инструментов KU 312 / KU 313, 500 200 100 мм, с ручкой / без ручки

ИМ-7.107461 до 22.04.2024

сталь, прямоугольный, для инструментов, длина 500 мм, с крышкой

Nora instruments Medizintechnik GmbH, Германия



Lid height		
16,5 x 8,0 x 3,5 cm	KU 31216	
20,0 x 10,0 x 5,0 cm	KU 31220	
22,0 x 15,0 x 5,0 cm	KU 31222	
16,5 x 15,0 x 5,0 cm	KU 31216B	
20,0 x 20,0 x 5,0 cm	KU 31220B	
22,0 x 15,0 x 5,0 cm	KU 31222B	
16,5 x 20,0 x 10,0 cm	KU 31216S	
20,0 x 20,0 x 10,0 cm	KU 31220S	
22,0 x 15,0 x 10,0 cm	KU 31222S	
Tray lid height		
16,5 x 8,0 x 3,5 cm	KU 31316	
20,0 x 10,0 x 5,0 cm	KU 31320	
22,0 x 15,0 x 5,0 cm	KU 31322	
16,5 x 15,0 x 5,0 cm	KU 31316B	
20,0 x 20,0 x 5,0 cm	KU 31320B	
22,0 x 15,0 x 5,0 cm	KU 31322B	
16,5 x 15,0 x 10,0 cm	KU 31316S	
20,0 x 20,0 x 10,0 cm	KU 31320S	
22,0 x 15,0 x 10,0 cm	KU 31322S	